


I ntraocular use of biologics

- The intraocular use of biologics has revolutionised the treatment of blinding diseases, especially in our ageing population.
- Preclinical development is hampered by anti-drug antibodies (ADAs) that form in animal models against human candidate medicines. ADAs limit development of products designed for longer duration of action.
- The PK-Eye™ replicates the human clearance times of intraocular protein medicines.
- The PK-Eye™ makes it more feasible to determine intraocular clearance times, dose escalation PK effects and stability profiles of biologics during preclinical development.

The PK-Eye™ replicates the human clearance times of biologics

- Biologics such as antibody based medicines and other proteins are large, charged molecules that are NOT permeable through the back of the eye (RCS pathway).
- Aqueous outflow is through the front of the eye (anterior chamber) and is the primary means of mass transfer in the eye needed to nourish the cornea and lens.
- Protein-based medicines primarily clear from the front of the eye (anterior pathway) due to aqueous outflow. Large molecules diffuse more slowly than small molecules in the vitreous.



- ❑ Using simulated vitreous (SV), the PK-Eye™ shows the clearance half-life ($T_{1/2}$) for ranibizumab (Lucentis®) is very similar to what is observed in humans.

Ranibizumab dose of 0.5 mg, 50 uL	
The PK-Eye™ (SV)	Humans
8.0 ± 3.1 days	7-9 days

